



# Gold Thin Films

## What is a gold thin film?

- A nanometer-scale film of gold deposited onto a supporting substrate

## What are gold thin films used for?

- Gold thin films have diverse applications in life science research, sensor development, reflectivity analysis, surface plasmon resonance detection, atomic force microscopy, and more

## What are the benefits of gold thin films?


- Chemically inert surfaces, resistant to oxidation
- Excellent conductors of electricity
- Readily functionalized, especially with alkanethiol monolayers
- Available with optional atomically smooth surface as “**ultra-flat**” gold




# Platypus Gold Thin Film Products



**Microscope slides**  
10 nm, 50 nm, 100 nm  
gold films



**Coverslips**  
22 mm square:  
10 nm, 50 nm gold films  
15 mm round: 10 nm gold film



**Silicon wafers**  
100 nm gold film




**Mica**  
200 nm gold film



**Ultra-Flat Films**  
100 nm on glass  
200 nm on mica



# Platypus Custom Gold Thin Films



**Custom gold and titanium coatings  
on your substrates of glass, silicon,  
or mica**

**Substrate**

Minimum 10 x 10 x 0.1 mm (0.25 x 0.25 x 0.0025 in)

**Dimensions**

Maximum 165 x 165 x 3 mm (6.5 x 6.5 x 0.12 in)

**Substrate Type**

Vacuum-compatible to  $10^{-7}$  Torr

**Gold Thickness**

Minimum 0 nm, Maximum 300 nm

(<7.5 nm may not provide a continuous coating)

**Titanium Thickness**

Minimum 0 nm, Maximum 50 nm


(<7.5 nm may not provide a continuous coating)



# Platypus® Ultra-Flat Gold Films

**Template Stripped Gold Chips** provide an extremely smooth and clean gold surface for a variety of research needs, including AFM, SEM, sensor development, and surface chemistry studies.

**AFM Characterization:** RMS roughness over a randomly selected 1  $\mu\text{m}^2$  area



1  $\mu\text{m}^2$  AFM scan of as-deposited gold on Silicon Wafer  
RMS roughness = 17.3 Å

1  $\mu\text{m}^2$  AFM scan of Platypus® Template Stripped Gold – Silicon Wafer  
RMS roughness = 3.6 Å

1  $\mu\text{m}^2$  AFM scan of Platypus® Template Stripped Gold – Mica  
RMS roughness = 4 Å

**Goniometer Characterization:** Water static contact angle

Surface Condition	Static Contact Angle
As-Deposited Gold (fresh)	25 - 40°
As-Deposited Gold (after exposure to atmosphere for 5 hours)	75 - 90°
Template Stripped Gold Chip (stripped prior to reading)	25 - 40°

*Note: All measurements carried out at ambient temperature.*



# Platypus Custom Glass Scribing

## Accurate, diverse shapes for flat substrates:


	Inches	Millimeters
Thickness	0.004 - 0.2	0.10 - 5.00
Max sheet size	24 x 24	600 x 600
Rectilinear accuracy	± 0.002	± 0.0508
Shape accuracy	± 0.003	± 0.0762
Minimum curve radius	0.1	2.5

For substrates < 1 x 3 inches that require gold coating, we recommend coating larger sheets and scribing to size after coating



# Gold Coating Methods

- **Thermal vapor deposition:** The material to be deposited is evaporated by electrical resistive heating at low pressure and is condensed onto the substrates
- **Sputtering:** The material to be deposited is ejected from a source by plasma discharge which is then deposited onto the substrates
- **Electron-beam vapor deposition:** The material to be deposited is heated to a high vapor pressure by electron beam bombardment in an ultra-high vacuum and is condensed onto the substrates. This method is used by Plastypus because it produces the smoothest, cleanest gold surfaces





# Advantages of Platypus Gold Films

- Ultra Clean Substrates** Substrates are cleaned by oxygen plasma prior to coating to promote adhesion and cleanliness
- Dedicated E-beam Coater** Only titanium and gold are used for coatings to ensure that the gold thin films are as pure as possible
- Cleanroom Facility** The gold coatings are applied in a class 10,000 cleanroom to reduce particulate contamination
- High Reproducibility** Dedicated gold coating equipment, a 99.999% pure gold source, calibrated deposition control, and consistent quality control measurements all contribute to the high reproducibility
- Titanium Adhesion Layer** Titanium is a better adhesion layer than Chromium. Chromium can diffuse into the gold film faster, contaminating it
- Low RMS Roughness** Roughness of Platypus' gold coated substrates will depend on the thickness of gold and type of substrate it is on, but generally the RMS roughness is between 1.5 and 4nm
- Custom Coatings** Platypus can work with customers to coat custom substrates with specified titanium and gold thicknesses



# Critical Parameters

## Thickness of gold layer

For most applications a 10 - 100nm thick gold thin film is used.

10-nm gold films are transparent whereas 100-nm gold films are opaque

## Roughness

A root mean square (RMS) roughness  $< 4$  nm is desired for most applications.

However, for certain scanning probe applications a lower roughness is required, which Platypus offers as **ultra-flat gold films**

## Surface purity

A clean surface is required for applications that employ gold surface modifications. Depending on the deposition process used, level of vacuum, purity of the gold source, and history of metals used in the coating equipment, the level of gold thin film purity can vary greatly between suppliers. Platypus uses only gold, titanium and vacuum-compatible substrate materials to ensure maximum cleanliness of gold thin films





# Cleaning Your Gold Films

## Atmospheric Deposition

Our gold thin films are fabricated in a clean room. However, it is challenging to prevent atmospheric components, especially organic compounds, from adsorbing to the highly active gold surface once they are in a typical laboratory atmosphere. We suggest users clean the gold immediately before use.

## Cleaning Gold Thin Films

If you plan to create self-assembled monolayers of alkanethiols on the gold, there is no need to clean the surface. The attachment of the alkane thiols is energetically favored over adsorbed atmospheric contaminants, and the alkanethiols will displace these.


Otherwise, a brief exposure to a butane flame will clean the surface of contaminants. We like **this butane torch model**, it is easy to use and after you have done your experiment you can use it to make crème brûlée.



# Glossary of Terms

**Nanometer (nm):**  $1 \times 10^{-9}$  m

**Ångström (Å):** 0.1 nm



**E-beam Evaporation:** A physical vapor deposition technique used for thin film coatings.

**Adhesion Layer:** A thin layer deposited for improving the adhesion of gold to the substrate.

**Roughness:** A measure of surface smoothness.

**Grain Size:** The size of gold grains present on the surface. Grain size impacts roughness.

**Self Assembled Monolayers (SAMs):** A single layer of molecules, most commonly alkanethiols, that spontaneously assemble into monolayers on a gold surface.



# QUESTIONS?

## **Platypus Technologies, LLC**

5520 Nobel Drive, Suite 100

Madison, WI 53711

T: 608.237.1270

Toll-Free: 1.866.296.4455

Fax: 608.237.1271

**[info@platypustech.com](mailto:info@platypustech.com)**



**[www.platypustech.com](http://www.platypustech.com)**